



# Gel Knit® Removable Dressings

Removable nasal and sinus dressings

## Riemann, Goodman, and Mannheim Gel Knit nasal and sinus dressings:

### Promotes optimal surgical outcomes

- Provides a moist physical barrier for tissue opposition and prevents the formation of adhesions post operatively

### Controls minor bleeding

- Moist and self-lubricating Carboxymethylcellulose – (CMC) infused fabric minimizes patient discomfort upon removal

### Allows for easy and early removal

- CMC-infused fabric remains lubricious for easy removal at 24-72 hours without adhering to nasal tissue



Soak dressing to hydrate Gel Knit hydrocolloid fabric

## Available in a range of sizes:



### Riemann Sinus Dressing *Rapid Rhino® 300 and 400*

For ESS procedures, Polypectomy

- Available in two sizes:
  - 3.0 cm (30 mm length x 6.5 mm width)
  - 4.0 cm (40 mm length x 6.5 mm width)
- 20 units per box



### Goodman Nasal Dressing *Rapid Rhino 500*

For Septoplasty, Turbectomy, Polypectomy

- Available in 5.5 cm (55 mm length x 5 mm width)
- 20 units per box



### Mannheim Nasal Dressing *Rapid Rhino 800*


For interior Septoplasty, Turbectomy, Polypectomy

- Available in 8.0 cm (72 mm length x 6.3 mm width)
- 20 units per box


# Gel Knit Removable Dressings

## Removable nasal and sinus dressings

### A randomized trial of Rapid Rhino Riemann packs following endoscopic sinus surgery<sup>1</sup>



### In a side-by-side comparison against Merocel®, patients reported significantly less pain on removal of the Gel Knit Nasal Dressing<sup>2</sup>



#### References

- Cruise, A.S., Amonoo-Kuofi, K., Srouji, I., Kanagalingam, J., Georgalas, C., Patel, N.N., Badia, L., & Lund, V.J. Clin. A randomized trial of Rapid Rhino Riemann packs following endoscopic sinus surgery. *Otolaryngology* 2006, 31, 25-32.
- Arya AK, Butt O, Nigam A. Double-blind randomized controlled trial comparing Merocel with Rapid Rhino nasal packs after routine nasal surgery. *Rhinology*, 2003; 41: 241-243.



#### International

**ArthroCare ENT**  
Skeppsbron 2, 111 30 Stockholm, Sweden  
Tel: +46 8 546 172 00  
Fax: +46 8 546 172 39  
ent@arthrocare.se

#### Australia

**ArthroCare (Australasia Pty Ltd)**  
PO Box 2186, Normanhurst NSW 2076  
Australia  
Tel: +61 2 9479 9200 Fax: +61 2 9479 9299  
service@arthrocare.net.au

#### Austria

**Atlantech Med. Prod. Vertriebs-AG**  
Kaufmannstrasse 15, 6020 Innsbruck, Austria  
Tel: +43 512 36 15 38-00 Fax: +43 512 36 15 38-60

#### France

**ArthroCare France**  
9, rue Noirot, 70000 Vesoul, France  
Tel: +33 84 76 81 38 Fax: +33 3 84 76 74 63

#### Germany

**ArthroCare (Deutschland) AG**  
Berghauser Str. 62, 42859 Remscheid, Germany  
Tel: +49 2191 93342 0 Fax: +49 2191 93342 49  
info@arthrocare.de

#### Scandinavia

**ArthroCare Scandinavia ApS**  
Science Center Skejby, Brendstrupgårdsvej 21  
DK-8200 Århus N., Denmark  
Tel: +45 7023 2850 Fax: +45 7023 2851  
info@arthrocare.dk

#### UK

**ArthroCare UK Ltd**  
St James Business Park  
9, Grimbold Crag Court, Knaresborough  
North Yorkshire HG5 8QB, United Kingdom  
Tel: +44 (0) 1423 888806 Fax: +44 (0) 1423 888816  
info@arthrocareuk.com

#### USA

**ArthroCare ENT**  
680 Vaqueros Avenue  
United States  
Sunnyvale, CA 94085-3523  
Tel: +1 408-736-0224

Not for distribution in the United States.

© 2009 ArthroCare Corporation. ArthroCare, Rapid Rhino, and Gel Knit are trademarks and/or registered trademarks of ArthroCare Corporation. Merocel is a registered trademark of Medtronic Xomed.

P/N 16717.gb Rev B