

VAC 30/45

The Mobile Suction Unit

Suitable for
disposable and reusable
liquid containers



otopront®

VAC 30/45

The VAC 30/45 is a very powerful unit, that is suitable for all kinds of suction procedures. Due to its easy mobility on antistatic rubber castors, the VAC 30/45 can be used in a variety of locations, the hospital clinic, on the ward for bedside visits, or indeed in the specialist's consulting room. The low noise and maintenance free pump is designed for continuous operation and reaches high performance values. The suction process either starts via photodetector or a footswitch underneath the OR-table. Beside the vacuum indicator and vacuum gauge, the unit is equipped with an autoclavable suction pipe cleaner, double overflow control fuse and bacterial air filter. According to the latest standards, the unit is available with disposable or reusable bottles, to collect the secretion. Tube connections are either snap systems or olive systems, that suit to all standard suction tubes and thus different kinds of currettes and canulas.

Technical Data

Suction capacity	95% max. vacuum, 35/45 liters /minute
Liquid container	disposable or reusable containers
Container capacity	disposable 1.000 to 3.000 ml reusable 1.500 ml, 2.500 ml, 5.000 ml
Level control	double overflow fuse
Tube connection	snap system or olives connection
Cleaning	suction pipe cleaner, autoclavable adapter
Cleaning reservoir	1. liter bottle for disinfection liquid
Suction tubes	different length and diameters available, automatic on/off switch via photodetector or footswitch
Casing	Steel casing with 2- layer paint
Mobility	4 antistatic rubber castors
Weight	28 kg
Dimensions	500 x 980 x 500 mm (WxHxD)
Power supply	110 V/230 V, 50/60 Hz
CE mark	accord. to MDD 93/42/EEC



The liquid- and the disinfection bottle are easy to change.



The double overflow fuse of the VAC allows a safe function of the silent vacuum pump.

All designs and specifications subject to change without notice.

Happersberger otopront GmbH
Langgasse 90
D-65329 Hohenstein

Tel. +49 (0) 61 20 - 92 17 0
Fax. +49 (0) 61 20 - 92 17 60

www.otopront.com
E-Mail: info@otopront.com